

# Brush head replacement

- 1. To remove:** Place QuickRelease™ tool\* on neck of toothbrush as shown and press (A), then slide brush head off toothbrush (B).
- 2. To replace:** Lay new replacement brush head parallel to neck of toothbrush (C), then slide head toward handle until lock clicks into place (D). Check that brush head is firmly attached. Do not insert brush head at angle (E).

Change your brush head when the blue Indicator® bristles fade halfway or every 3 months, as recommended by dentists or dental hygienists.

\* included in CrossAction® Power replacement brush.

